


StandUp Walker



StandUp Walker

- 2-in-1 Stand-Assist + Walker
- Weight capacity: 400 lbs
- 18" and 22" width options
- Height adjustable



Made in USA



URise Products, LLC
Charlotte, NC

Designed and Assembled in USA

www.URiseProducts.com

1-844-99-URise

To learn more about StandUp Walker, including warranty information, visit our website at www.StandUpWalker.com

Please read all instructions before use.

StandUp Walker.com

18" & 22" User Guide



This is a user guide for both the 18" width and 22" width products. The 2-in-1 Stand-Assist + Walker is designed primarily for indoor use.

PRODUCT SPECIFICATIONS

Weight Capacity	400 lbs.
Product Weight	10.5 lbs.
Size Options	18" W and 22" W
Frame Warranty	Limited Lifetime

Get Up And Go With Confidence


URise Products has created video tutorials that match the descriptions below at www.StandUpWalker.com/use.

How to Fit the StandUp Walker

Step 1. When the user is standing, adjust the height of the handles so that they are even with wrists of the user.

Step 2. To adjust the handle height, first unlock the snap collars, then press the spring-loaded push button in and adjust the height up or down.

Step 3. Once the arm has been properly adjusted, be sure to secure the snap collar for stability.



How to Fold the StandUp Walker

Step 1. Release the snap collar on each arm. Press the corresponding push button and turn the telescoping arm so that the handle is towards the middle of the frame. Then secure the snap collar. This will keep the arms from rotating during transport. Repeat with second arm.

Step 2. To fold the legs, turn the StandUp Walker sideways on its frame. Push the button on the leg and rotate the leg 180 degrees to the corresponding hole location. Repeat with other leg. Each leg should contour the wheel for a compact travel fit.

PRECAUTIONS:

Do not use on stairs, moving sidewalks or escalators.

BE AWARE OF OBJECTS AND OBSTRUCTIONS DURING USE.

How to Use the StandUp Walker

While in a seated position, be sure the StandUp Walker is within reach of the user.

Step 1. To engage for use, lift up on each of the two levers in the middle of the unit simultaneously to release the security pins. The StandUp Walker will then drop to the lower ergonomic position.



Step 2. Pull the walker as close as possible so that the skis are positioned under the seat. There are three different hand positions: pistol, medium and the recommended ball-joint grips.



Step 3. As the user rises, their weight can be placed entirely on the StandUp Walker.



Step 4: Once standing, lift the StandUp Walker using the ball-joint grips attached to the handles. The hinges will automatically secure when lifted into the standing position. Listen for the audible click sound before engaging the StandUp Walker for mobility.

