

Clever-Lite Walker, Adult, with 8" Casters

HCPCS: E0143+E0156



FIGURE A

FIGURE B

- Front wheels can be set in either "swivel" (Figure A) or "fixed" (Figure B) position
- Special loop lock made of internal aluminum casting operates easily and ensures safety
- Sturdy 1" diameter anodized aluminum construction
- Comes with flip-up seat with built-in carry handle
- Allows individual to be seated or use as traditional walker
- Folds with patented side-paddle release

PRODUCTS	Description
10230BRAKE	HandBrake&CableAssembly,10230,1/ea
10243	Clever-Lite Rollator,w/8 wheel,1/cs

SPECIFICATIONS	Description
Carton Shipping Weight	21 lbs
Depth (opened at base)	25.5"
Handle Height	35" - 39"
Seat Dimensions	19" (W) x 11.25" (D)
Seat Height	24" - 29"
Warranty	Limited Lifetime
Weight	18 lbs
Weight Capacity	300 lbs
Width	18.5" inside hand grip 21" inside back legs 25.5" Overall

PARTS	Description
10243F	Front Fork w/tubing,10243,1/ea, EA
10230SF	Seat Frame,10230,1ea, EA
10243RW	Rear Wheel for 10243,1/ea, EA
10243W	Front Wheel for 10243, 1/ea, EA
10230BRAKE	HandBrake&CableAssembly,10230,1/ea, EA
12002PB	Push Button for Bath Bench,Walker,1/ea, EA
9500001232	Brake Cable/Adjust,Lower,1/ea, EA
10230SN	Blow Molded Seat,10230,1ea, EA
9500075010	Front Axle Bolt,750,10243,1/ea, EA
10230B	Basket,Clever-Lite,Adult,1/ea, RTL, EA
950340129901	Cable,10230,1/ea, EA