

Deluxe 3 Wheel Aluminum Rollator, 7.5" Casters

HCPCS: E0143



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F

- Comes standard with basket, tray and pouch (Figure A)
- Basket conveniently folds for easy storage (Figure B)
- Ships in ultra-compact carton Comes with frame folded and compact Pull lightweight, aluminum trombone-style frame to lock in place and secure with lock knob
- 7.5" casters for indoor or outdoor use
- Handles easily adjust in height with self-threading knob (Figure C)
- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure D)
- Easy, one-hand folding

PRODUCTS

Description

10289BL	3 Whl Alum Rollator,Blue,1/ea
10289RD	3 Whl Alum Rollator,Red,1/ea

SPECIFICATIONS

Description

Carton Shipping Weight	17 lbs.
Handle Height	31" - 38"
Unit Dimensions	25" (W) x 24" (L)
Warranty	Limited Lifetime
Weight	11 lbs.
Weight Capacity	300 lbs

PARTS

Description

9503E17101-L	Hand Brake with Cable, 171, Left, 1/ea, EA
9503E17101-R	Hand Brake with Cable, 171, Right, 1/ea, EA
816T	Tray, RTL199, 1/EA, EA
9503E1028901-B	Handlebar,BlackFinish,10289,Right/Lft,1e, EA
815B	Folding Basket, RTL199, 1/EA, EA
808	BLACK TOTE,1/ea, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
9503E17109	Fork,171,1ea, EA
9503E17106	Rear Axle with Nut,171,1ea, EA
9503E17107	Front Axle with Nut,171,1ea, EA
9503E17108	Wheel Assembly,171,1ea, EA
9503E17111	Upper Pivot,171,1ea, EA
9503E17112	Folding Hinge,171,1ea, EA
9503E17110	Lower Pivot Bolt,171,1ea, EA

PARTS

Description

9503E17103	Height Adjustment Knobs, RTL199,1/PR, PR
9503E17104	Wheel Brake,Left,171,1ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
950330049901	Hand Grip, Right, 3-wheel Rollator, 1/ea, EA
950330039901	Hand Grip, Left, 3-wheel Rollator, 1/ea, EA