

Durable 4 Wheel Rollator with 7.5" Casters

HCPCS: E0143+E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Comes with durable, padded seat
- Removable, hinged, padded backrest can be folded up and down
- Comes with a basket to store personal items
- Durable steel frame
- 7.5" Non-marring casters are ideal for indoor and outdoor use
- Easy-to-use deluxe loop locks (Figure A)
- Lever allows for independent handle height and angle adjustment (Figure B)

PRODUCTS	Description
10257BL-1	Rollator,Steel,4 Wheel,Blue,1/ea
10257RD-1	Rollator,Steel,4 Wheel,Red,1/ea

SPECIFICATIONS	Description
Carton Shipping Weight	22.2 lbs.
Handle Height	31" - 37"
Seat Dimensions	12" (L) x 12" (W)
Unit Dimensions	23.5" (W) x 25.5" (L)
Warranty	Limited Lifetime
Weight	18.6 lbs.
Weight Capacity	300 lbs

PARTS	Description
9502F1025714	Right Handlebar, 10257, 1/ea, EA
9502F1025715	Left Handlebar, 10257, 1/ea, EA
95032F260N06	Dust Cap,260N,10257, 1/pr, PR
9502F1025712	Padded Seat W/Brackets, 10257, 1/ea, EA
9502F1025702	SeatBracketsW/Screws,HardSeat,10257,1/pr, PR
9502F1025713	SeatBrackets&Hrdwr,PaddedSeat,10257,1/pk, PKG
9502F1025711	WheelBrakeRight10257,1ea, EA
9502F1025710	WheelBrakeLeft10257,1ea, EA
9502F1025707	FoldingHingeAssembly10257,1ea, EA
9502F1025709	FrontAxleBolt10257,1ea, EA
95032F260N07	Stem Bearing Set,260N,1Set EA, EA
95032F260N14	Right Hand Grip,260N,1ea, EA
95032F260N15	Left Hand Grip260N,1ea, EA
95032F260N10	Height Adjustment Knobs,260N, 10257 1pr, PR

PARTS**Description**

95032F260N09	Wheel W/Bearings, Grey, 260N & 10257,1ea, EA
95032F260N09-BLK	Wheel W/Bearings, Black,260N & 10257,1ea, EA
9502F1025708	Fork10257,Blue,1ea, EA
9502F1025708R	RedFork10257,1ea, EA
95032F260N12	Right Handlebar,260N,1ea, EA
95032F260N13	Left Handlebar,260N,1ea, EA