

Aluminum Rollator, 7.5" Casters

HCPCS: E0143+E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G




FIGURE H

- Comes with new seamless padded seat (Figure A)
- Zippered pouch under seat for added privacy and security of personal items (Figure B)
- 7.5" Black non-marring tires are ideal for indoor and outdoor use
- Removable, hinged, padded backrest can be folded up or down
- Easy-to-use deluxe loop locks (Figure C)
- Brakes with serrated edges provide firm hold (Figure D)
- Ergonomic handles are easy-to-grip and are adjustable

PRODUCTS

Description

R728BL	Rollator,7.5" 4Whl,Alum,Blue,EA
R728RD	Rollator,7.5" 4Whl,Alum,Red,EA

SPECIFICATIONS

Description

Carton Shipping Weight	17 lbs.
Handle Height	33" - 38"
Seat Dimensions	14" (W) x 23" (H) x 12" (D)
Unit Dimensions	24" (W) x 28" (L)
Warranty	Limited Lifetime
Weight	14 lbs.
Weight Capacity	300 lbs

PARTS

Description

9505W1026120	Tote,10261,726,728,R800,R900,1/ea, EA
10208KN	Knob,Ht.Adj.Rollator,w/Attch.1/ea, EA
9501S72605R	FullBrakeSysHndGrips&HndTubeRT,726,1/ea, EA
9501S72605L	FullBrakeSysHndGrips&HndTubeLFT,726,1/ea, EA
9501S72604	SeatW/Brackets&AttachHardwre,726/728,1ea, EA
9501S726BK	Back Rest with Foam,Black,726,728, 1/ea, EA
9500001232	Brake Cable/Adjust,Lower,1/ea, EA
750CLIP	Wheel LockClip,750N,1/ea, EA
9500004101	Cable,Fastner,1/ea, EA
9501S72800-01	Front Wheel for 728/R728,1/ea, EA
950340059901-2M	Padded Seat Brackets, 2M Ser.#, 1/pr, PR
9502M72801	FrontFork/CasterAssembly,728,2M Ser.#,1ea, EA
9501S72603	Backrest Pivot Brackets for R726,1 pr, PR
9501S72607	Brake System Only for 726,1/ea, EA

PARTS	Description
9501S72616	Spring,R726,1/ea, EA
9500075011	Rear Axle Bolt,750,780,1/ea, EA
9501S7261L	Wheel Lock, Left, R726,1/ea, EA
9501S7261R	Wheel Lock, Right, R726,1/ea, EA
9501S72800	Rear Wheel for 728/R728,1/ea, EA
STDS1002-1	Axle for 8" caster, 1/ea, EA