



EXCELTRANSLATOR™

The most comfortable combination Rollator/Wheelchair on the market!



Walk,



Rest,



Ride!

New!
Patent Pending

Combine the functionality of both a rollator and wheelchair into one unit!

In just seconds, our Excel Translator easily converts from a rollator to a wheelchair and vice-versa. The Excel Translator allows rollator users to sit down and ride in a wheelchair when their legs get tired. And, wheelchair users can stand up and walk along when they want to get more exercise. Best of all, you only need to own one piece of equipment, rather than two!

The Excel Translator weighs only 18 lbs which is no heavier than the average rollator or transport chair. It comes with a convenient cup holder, under-the-seat basket, zippered side carrying case, and easy-to-adjust push-button foot rests. You'll be delighted with the comfortable, spacious seat which is more than 75% roomier than our average rollator seat.

Excel Translator – The most comfortable Combination Rollator/Wheelchair on the market!



- Height-Adjustable Push Grips

- Breathable Nylon Upholstery

- Comfortable Hand Brake

Patent Pending!

- Easy-to-Adjust Pushbutton Footrests Lock to Sides of Translator When Not In Use

Flip it to convert it!



- Convenient Side Carrying Case
- Convenient Cup Holder
- Roomy Under-Seat Basket
- Strong Carrying Handles



Specifications:

Item Number: MDS808200TR

Seat Size: 19" (48.3 cm) W
x 16" (40.6 cm) D

Seat Height: 19" (48.3 cm)

Frame Color: Cobalt Blue

Weight Capacity: 250 lbs (113.4 kg)

Unit Weight: 18 lbs (8.2 kg)

(does not include footrests or cases)

Casters: 6" (15.2 cm) front

7" (17.8 cm) rear

HCPCS: Pending



Medline Industries, Inc.

One Medline Place
Mundelein, IL 60060

Medline United States

1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada

1-800-396-6996

www.medline.ca | canada@medline.com

Medline México

01-800-831-0898

www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. Please contact your Account Representative for additional information.

We reserve the right to correct any errors that may occur within this flyer. ©2009 Medline Industries, Inc. Medline is registered trademark of Medline Industries, Inc.

Excel Translator™ is a registered trademark of Medline Industries, Inc.

MKT209192 / LIT762 / 10M