

Mobility

Free2Go Rollator

The Free2Go Rollator combines the benefits of the rolling walker (“rollator”) with a toilet safety product, thereby allowing one to safely, confidently and discreetly use a toilet in **any** restroom.

Upon entering a restroom, push the rollator over the toileting area, engage the wheel locks and raise the padded seat to use. The flip-up padded seat is equipped with a discreet pouch to conveniently store cleaning wipes for after each use.



Walking Assist Mode

Toilet Assist Mode

FEATURES

- Combines the benefits of a rollator, toilet safety frame and raised toilet seat
- The frame of the rollator provides support when sitting and standing
- Provides additional height over a standard toilet to reduce body strain
- Restores dignity and independence
- Reduces the fear of falling and the possibility of life threatening accidents
- Portable, lightweight and discreet design
- 8" flat wheels for added stability
- Removable flip-up curved padded backrest
- Vinyl padded seat with easy-folding system
- Storage pouch located under seat ideal for cleaning wipes
- Limited lifetime warranty on frame and 90 days on parts

Product Specifications:

Seat Measurements:	15" (W) x 15 ½" (D)
Toilet Seat Opening:	8 ¼" x 9 ½"
Seat-to-Floor Height:	23 ½" (top of seat) 22 ¼" (top of toilet seat)
Adjustable Handle Height:	33 ½" – 38"
Unit Width:	24"
Product Weight:	20.2 lbs.
Patient Weight Capacity:	300 lbs.

Item #	Description	Color	Qty
F2G-ROL8CBG	Free2Go Rollator with Toilet Seat, 8-inch Wheels	Burgundy	Each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com