

TRANSITIONS®



EZ-ACCESS®

Threshold ramps, mats, and plates

EZ-ACCESS® knows that a gentle slope can make a world of difference. TRANSITIONS® threshold solutions allow for a smooth ground-to-sill transition over doorway thresholds throughout your home. Reduce preventable trip hazards and make your home more comfortable and safe for mobility users.



The TRANSITIONS® Advantage


- Sturdy, seam-free surface makes TRANSITIONS ramps, mats, and plates the safe solution for wheelchair, scooter, walker, and cane users
- Attractive designs are available in multiple styles and materials to suit your needs and taste
- Slip-resistant tread and surfaces
- Durably designed for years of continued use
- Quick and easy setup and removal - installs in minutes so there's little-to-no disruption for you and your family
- Ideal for doorways, sliding glass doors, and raised landings
- Portable, convenient, and durable
- The entire line is backed by our lifetime warranty
- Made in the USA




Visit ezaccess.com for product specifications

TRANSITIONS® Modular Entry Ramp

- Modular design with interchangeable extensions can accommodate thresholds from 3/4 to 6-1/4 inches high
- Installation is easy with our hidden hinge ... simply add and remove extensions to accommodate your needs
- Pivoting flap attaches to the top of the ramp to provide 1/4" of additional adjustability
- Made of aluminum and features an extruded, slip-resistant surface for superior traction
- 34 inches wide
- 700-pound weight capacity



TRANSITIONS® Angled Entry Ramp



- Swivel feet accommodate varying heights, angles, and surface types
- Three lengths are available (12, 24, and 36 inches) to accommodate thresholds ranging from 1-3/8 to 5-7/8 inches high
- Manufactured using high-strength aluminum that will not rust or corrode
- Features an extruded, slip-resistant surface for superior traction
- 36 inches wide
- 700-pound weight capacity

TRANSITIONS® Angled Entry Mat



- This "welcome" mat - made using recycled rubber - is eco-friendly and can be used safely on nearly any surface
- Stands up well in harsh outdoor conditions
- Features beveled edges for cross-traffic access
- Available in 1-1/2 and 2-1/2 inch heights and can easily be trimmed or notched if needed
- Usable width of 36 inches (40-43 inches overall)
- 850-pound weight capacity
- Available in three colors: Black, Storm Grey, and Hazelnut Brown

TRANSITIONS® Angled Entry Plate

- Formed aluminum construction for simple, one-step threshold navigation
- Three lengths are available (8, 10, and 12 inches) to accommodate thresholds up to 2 inches high
- Features an applied, slip-resistant surface for greater traction, even when wet
- 32 inches wide
- 700-pound weight capacity






TRANSITIONS® Angled Entry Ramp

Technical Specifications

MODEL	TAER 12	TAER 24	TAER 36
OVERALL DIMENSIONS	11.877" x 36.25"	23.875" x 36.25"	35.904" X 36.25"
HEIGHT (RAMP ONLY)	1.362" (1-3/8")	1.669" (1-11/16")	3.282" (3-5/16")
USABLE HEIGHT W/ 1" SWIVEL FEET	1.92" TO 2.958" (1-15/16" TO 2-3/4")	N/A	N/A
USABLE HEIGHT W/ 1½" SWIVEL FEET	N/A	2.154" TO 4.361" (2-1/8" TO 3-3/8")	N/A
USABLE HEIGHT W/ 2½" SWIVEL FEET	N/A	3.178" TO 4.361" (3-3/16" TO 4-3/8")	3.681" TO 5.836" (3-5/8" TO 5-7/8")
RATED LOAD (DO NOT EXCEED)	700 lbs., 350 lbs. per caster or wheel	700 lbs., 350 lbs. per caster or wheel	700 lbs., 350 lbs. per caster or wheel
WEIGHT	5.15 lbs.	10.4 lbs.	15.15 lbs.
WARRANTY	Lifetime Limited	Lifetime Limited	Lifetime Limited


Weights and dimensions are approximate and subject to change.

TAER 12



Shown with removable 1" swivel foot installed

TAER 24



Shown with removable 1½" swivel foot installed

TAER 36



Shown with removable 2½" swivel foot installed

LIFETIME WARRANTY. Please register at www.ezaccess.com/warranty-satisfaction.
© EZ-ACCESS®, a division of Homecare Products, Inc. All rights reserved.

All text and/or images contained in this document are proprietary and may not be shared, modified, distributed, reproduced, or reused without the express written permission of EZ-ACCESS.