

Seat2Go Positioning Seat

HCPCS: E1399-1,E1399



FIGURE A

- Contoured seat prevents child from sliding forward and improves seating posture
- 95° seat to back angle
- Double-strap trunk support and a positioning belt
- Two sets of straps secure Seat2Go to any chair
- When not in use, the seat straps can be stored under the seat behind the clear flap and the back straps in the back pocket
- Can be used as an insert for a single or double stroller
- Weighs only 5 lbs

PRODUCTS	Description
CRS 2000	Seat2Go, small,1/cs
CRS 3000	Seat2Go, Medium,1/cs

SPECIFICATIONS	Description
Carton Shipping Weight	8 lbs 11 lbs
Warranty	Limited 2 Year

ACCESSORIES	Description
CRS 8025	Abductor, Seat2GO sm/lg, 1/bx, BX
CRS 8000	Headrest, Seat2Go & Back2Go, 1/ea, BX