

MSS Tilt & Recline Seating System

HCPCS: E1399-1



FIGURE A



FIGURE B




FIGURE C



FIGURE D

- Multi-Axis Headrest - height, depth and angle adjustable
- It can be placed off center
- Seat - depth adjustable
- Back - height and angle adjustable
- Harness - adjustable zippered "H" harness for anterior chest support
- Pelvic Belt - adjustable in length
- Laterals - height and width adjustable
- They assist in stabilizing and aligning trunk and pelvis
- Each lateral is individually adjustable to control unwanted weight shifts
- Armrest - height and angle adjustable
- Adductors - two sizes are included to facilitate a range of chair widths
- The adductors control lower extremity positions and prevent adduction of user's legs
- Abductor - height and depth adjustable
- Footrests - height, depth and angle adjustable They include foot sandals with toe loop
- Tray - attaches onto armrest and is height adjustable
- It provides trunk stabilization and an easy-to-clean surface with a safety lip
- It is ideal for school, play and feeding
- Each footrest is individually adjustable for custom fit
- Removable and washable fabric covers

PRODUCTS

Description

| | |
|----------|---------------------------|
| MS 3000N | MSS Tilt and Recline,1/ea |
|----------|---------------------------|

SPECIFICATIONS

Description

| | |
|-------------|-----------|
| Back Height | 14" - 19" |
|-------------|-----------|

| | |
|------------------------|--------|
| Carton Shipping Weight | 40 lbs |
|------------------------|--------|

| | |
|---|----------|
| Distance Between Laterals (FC 2000N FC 4000N) | 9" - 15" |
|---|----------|

| | |
|----------------------|----------|
| Headrest Height Adds | 4" - 10" |
|----------------------|----------|

| | |
|---------------|-----|
| Overall Width | 18" |
|---------------|-----|

| | |
|------------|------------|
| Seat Depth | 3.5" - 12" |
|------------|------------|

| | |
|------------|------------|
| Seat Width | 9" - 10.5" |
|------------|------------|

| | |
|-------------------------|-----|
| Seat-to-Footrest Height | 14" |
|-------------------------|-----|

| | |
|----------|------------------|
| Warranty | Limited Lifetime |
|----------|------------------|

| | |
|--------|--------|
| Weight | 35 lbs |
|--------|--------|

| | |
|-----------------|---------|
| Weight Capacity | 100 lbs |
|-----------------|---------|

| PARTS | Description |
|---------|--|
| FC 2024 | Dining Tray,FC2000,1/bx, BX |
| MS 4011 | AbductorWithoutCover,MS 3000N, 1/ea, EA |
| MS 4012 | Abductor Cover, MS 3000N, 1/ea, EA |
| MS 4013 | ArmrestAssemblyLeft, MS 3000N, 1/ea, EA |
| MS 4014 | ArmrestAssemblyRight, MS 3000N, 1/ea, EA |
| MS 4015 | ArmPad&Hardware, MS 3000N, 1/ea, EA |
| MS 4016 | Armrest Tilt Pin, MS 3000N, 1/ea, EA |
| MS 4017 | ArmrestBracket&Hardware, MS 3000N, 1/ea, EA |
| MS 4018 | L Shaped Footrest Bracket,MS 3000N, 1/e, EA |
| MS 4019 | FootrestMountingBracket, MS 3000N, 1/ea, EA |
| MS 4020 | Headrest Cover, MS 3000N, 1/ea, EA |
| MS 4021 | Back Cushion Only, MS 3000N, 1/ea, EA |
| MS 4045 | CoverforBackrestCushion, MS 3000N, 1/ea, EA |
| MS 4023 | Lateral w/o Cover, Right, MS 3000N, 1/e, EA |
| MS 4024 | Lateral w/o Cover, Left, MS 3000N, 1/ea, EA |
| MS 4025 | Lateral Pad Cover, MS 3000N, 1/pr, PR |
| MS 4027 | Footrest Rod, Right, MS 3000N, 1/ea, EA |
| MS 4026 | Footrest Rod, Left, MS 3000N, 1/ea, EA |
| MS 4029 | Foot Strap & Hardware, MS 3000N, 1/pr, PR |
| MS 4028 | FootSandal w/o MountBracket,MS 3000N,1pr, PR |