



Ultima Bath Chair
shown with optional head pad



FEATURES

- Independently adjustable seat, back and leg angles
- Seat-to-back angles: 0°, 30°, 60°, 90° (from vertical)
- Seat-to-leg angles: 0°, -30°, -60° (from horizontal)
- Bacteria and fungal resistant ComfortMesh™ seat cover
- Includes two positioning belts
- Folds flat for storage



Beach Bubble Blue



Lionfish Lavender



Ultima Access Bath Chair
Seat option provides easy access
for personal care

ACCESSORIES



Detachable Extension Base



Detachable Rolling Base



Attached Rolling Base*



Foldable Transfer Base*



Compact Transfer Base*



Head Support



Lateral Support Pads



Rapid-Dry Pelvic Positioning Belt



Pelvic Positioning Belt



Detachable Armrest

*See Bath Transfer literature for more information on Transfer Base Options



ORDER FORM FAX TO 310-305-1718 OR SUBMIT ONLINE

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

BASE UNIT	ULTIMA			ULTIMA ACCESS	
Beach Bubble Blue 	BCU-8620M-BB	BCU-8720L-BB	BCU-8920XL-BB	BCUA-8726L-BB	BCUA-8926XL-BB
Lionfish Lavender 	BCU-8620M-LL	BCU-8720L-LL	BCU-8920XL-LL	N/A	N/A
	\$984.50	\$984.50	\$984.50	\$1391.50	\$1441
Product Size	Small/Medium	Large	Extra Large	Large	Extra Large
Weight Capacity	130 lb	160 lb	180 lb	160 lb	180 lb
Height Capacity	40" -54"	54" - 68"	68" - 80"	54" - 68"	68" - 80"
Overall Length (Flat)	45"	53"	61"	53"	61"
Chair Leg Length (Flat)	39.5"	41.5"	44.5"	41.5"	44.5"
Overall Width	17"	18"	20"	18"	20"
Inside Seat Width	14"	15"	17"	14"	16"
Seat Depth	11"	13"	15"	14"	17"
Calf Length	10"	13"	15"	13.5"	15.5"
Back Height	22"	27"	31"	26"	30"
Seat Height Above Floor	4" - 15"	4" - 15"	4" - 15"	4" - 15"	4" - 15"
Product Weight	20.5 lb	21.5 lb	23 lb	28 lb	29.5 lb

ACCESSORIES (fit all sizes of bath chair unless noted)

Detachable Extension Base (Raise chair by 9.5") 18"W x 24"-37"L x 9"H	\$770	BCUZ-A8880	
Detachable Rolling Base (Raise chair by 13") 25"W x 30"-44"L x 13"H	\$874.50	BCUZ-A8890	
Head Support Pads	\$170.50	BCZ-A8875-BB Beach Bubble Blue	BCZ-A8875-LL Lionfish Lavender
Lateral Support Pads	\$170.50	BCZ-A8830-BB Beach Bubble Blue	BCZ-A8830-LL Lionfish Lavender
Detachable Armrests	\$495	BTZ-A5051	
Positioning Belt	\$71.50	BCZ-A8835-G	
Pelvic Positioning Belt	\$93.50	BCZ-A8834	
Rapid Dry Positioning Belt (can be used as pelvic and/or lap belt)	\$269.50	BCZ-A8836M Fit range: 22" - 68"	BCZ-A8836L Fit range: 25" - 72"