

### 3-Position Recliner

HCPCS: E1031



FIGURE A



FIGURE B

- 3-positions: upright (Figure A), deep recline and elevated leg rest/footrest (Figure B)
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out
- Side panels "pop off" for easy cleaning and maintenance
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters (2 with locks) make it easy-to-move and secure the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Neutral-colored, sturdy steel frame is easy-to-clean and maintain
- Moisture barrier on seat prevents seepage
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117
- Standard

PRODUCTS	Description
D574-BR	3 Position Recliner,Blue Ridge,1/cs
D574-J	3 Position Recliner,Jade,1/cs
D574-R	3 Position Recliner,Rosewood,1/cs

SPECIFICATIONS	Description
Carton Shipping Weight	110 lbs.
Overall Dimensions	47" (H) X 25" (W) X 41" (D)
Seat Dimensions	21" (H) x 19" (W) x 19" (D)
Seat to Floor Arm Height	9" - 26"
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts
Weight Capacity	250 lbs
Width	20" between arms

PARTS	Description
D574P-1058BR	Footrest,D574,Blue Ridge,1/bx, BX
D574P-1058J	Footrest,D574,Jade,1/bx, BX
D574P-1055BR	Legrest,D574,BlueRidge,1/bx, BX
D574P-1058R	Footrest,D574,Rosewood,1/ea, EA
D574P-1055R	Legrest,D574,Rosewood,1/bx, BX
D574P-1055J	Legrest,D574,Jade,1/bx, BX
D574SEAT-BR	Seat Only,D574-BR,BlueRidge,1, BX
D574SEAT-R	Seat Only,D574-R,Rosewood,1/bx, BX

PARTS	Description
D574P-SEATLEG-J	Seat with LegRest,D574,Jade,1 ea, EA
D574P-SEATLEG-BR	Seat with LegRest,D574,Blue,1 ea, EA
D574P-SEATLEG-R	Seat with LegRest,D574,Rose,1 ea, EA
D574BACK-J	BackOnly,D574-J,Jade,1/bx, BX
D574BACK-BR	BackOnly,D574-BR,Bl.Ridge,1/bx, BX
D574BACK-R	BackOnly,D574-R,Rosewod,1/bx, BX
D574P-1048	HeadrestAdjBracket,D574/D574EW,Left,1 ea, EA
D574-BBN	Back Release Button,D574,1/ea, EA
D574P-1003	Hardware Pack,D574,1set/bx, BX
D574P-1061	Bolt Set,Footrest Sup.Bar,1/ea, EA
D574P-1036	Table Fitting,D574,1/ea, EA
D574P-1035	Recliner Mechanism,D574,Rt.&Lf,1/bx, BX