





The **Invacare® Semi-Electric Bed** combines effortless positioning of the upper body and knees with the economy of manual bed height adjustment. The pendant control provides motorized positioning of the upper body and/or knees. A manual crank located at the foot end of the bed adjusts bed frame height for caregiver assistance.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit **www.invacare.com/bedsafety** for additional safety and entrapment information.

Features

High impact bed end panels are durable, impact and scratch resistant; washable for easy cleaning

Universal bed end can be used as either the headboard or footboard on the bed

Rubber caster tread for better tracking with non-marring material

DC low voltage, compact motor system Quiet motor operation and fast cycle time on lowering height of bed for efficient bed set-up Motor system with waterproof enclosures is rated IPX6. Washable with standard water pressure. The bed must not be power washed

Color-coded connections on junction box for easier maintenance

Double-insulated junction box helps provide extra protection against electrical shock Ergonomic hand pendant comes with durable strain relief and heavy-duty cord. Waterproof casing for easy cleaning, conveniently attaches to bed sheet with clip or rail with hand pendant cradle (sold separately)

Hand pendant cradle available for convenient storage of hand pendant (model # 26975)

Interchangeability simplifies storage, reduces inventory complexity

HCPCS Code: E0295



- 1. Patented universal bed ends
- 2. Ergonomic hand pendant
- 3. Color coded connections
- 4. Emergency manual crank
- 5. Thermal overload protection
- 6. Double insulated electrical components/wiring

Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
5310IVC	88" L 15"-23" H 36" W	36" W	afe Working Loa 450 lb.* ax. Patient Weig 350 lb.		N/A	UL962 CSA Certified	5 Years Welds 2 Years Mechanical/ Electrical
5000IVC Headspring				45 lb.	43" L 8" W 36" D		5 Years Welds 2 Years Mechanical/ Electrical
5890IVC Footspring				77 lb.	47" L 9" W 36" D		5 Years Welds 2 Years Mechanical/ Electrical
5301IVC Universal Bed End				55 lb.	39" L 7" W		



<u>All references to HCPCS codes</u>: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

24" D

Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

©2017 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ^m and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 10-124 Rev. 03/17 170247

