

# NEPTUNE LIFT CHAIR

## Neptune CH4016 (Black), CH4017 (Brown) & CH4018 (Bone.)

The Neptune is a luxurious styled infinite position lift chair with very high quality genuine leather upholstery, studded arms and button back rest. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease movement of the chair.

With a Wide/Tall size rating, this chair will accommodate people of wider widths and above average heights.



Studded Arm Accents



Top Grain Cow Hide Leather



Button Back Design



### FEATURES:

- Lounger style that fully supports the legs (no gap between footrest and seat)
- Removable back for easy handling
- Easy to use hand control
- Maximum seat elevation with superb stability in full lift position
- Wide construction provides greater side to side stability
- Top steel base for durability and cross supports
- Engineered bushings and washers at pivot points provide quiet smooth operation and durability
- Direct drive activation synchronizes right and left sides for easy opening
- Low voltage motor, 24 volt, UL approved
- Battery backup
- Infinite position recline
- 4 small casters
- Genuine Leather Upholstery
- OKIN dual motor
- Foam meets CA117
- Button back



MODEL NAMES	NEPTUNE
Model numbers:	CH4016 (Black) CH4017 (Brown) CH4018 (Bone.)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	21.6" (55cm)
Seat Depth:	21.6" (55cm)
Back Height:	28.3" (72cm)
Seat Height:	19.7" (50cm)
Overall Width: (all positions)	35" (88.9 cm)
Overall Depth: (neutral position)	36" (91.4 cm)
Overall Height: (neutral position)	43" (109.2 cm)
Overall Depth: (reclined)	73.25" (186.1)
Overall Height: (reclined)	25.5" (68.4 cm)
Box 1	87cm x 90 x 64cm
Box 2	80cm x 80cm x 32cm
Net weight	122.45 lb (55.7 kgs) (box1) 28.65 lb (13.2 kgs) (box 2)
Gross Weight	126.45 lb (57.5 kgs) (box1) 30.65 lb (13.9 kgs) (box 2)