

PLUTO LIFT CHAIR



Pluto CH4007 (Blue) CH4008 (Brown)

The Pluto lift chair is an infinite position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair.

With a Medium/Average size rating, this chair will accommodate people of medium width and average height.



- Infinite position recline
- Four 2" casters included (rear casters lock)



The spring assisted footrest on all models will flip up if it comes in contact with an obstruction when lowering to minimize injury or damage.



Supplied locking swivel casters can be installed on all models, except Earth and Mars which have 3" casters. Ideal for users who may want to move their chair occasionally or rental units.

FEATURES:

- Lounger style that fully supports the legs (no gap between footrest & seat)
- Removable back for easy handling
- Easy to use hand control
- UL / CUL Approved
- Maximum seat elevation with superb stability in full lift position
- Wide construction provides greater side to side stability
- Top steel base for durability & cross supports
- Engineered bushings & washers at pivot points provide quiet smooth operation & durability
- Direct drive activation synchronizes right & left sides for easy opening
- Low voltage motor, 24 volt
- Battery Back Up



MODEL NAMES	PLUTO
Model numbers:	CH4007 (Blue) CH4008 (Brown)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	19.5" (49.5 cm)
Seat Depth:	19.5" (49.5 cm)
Back Height:	27.3" (69.3 cm)
Seat Height:	20.3" (51.6 cm)
Overall Width:	33.5" (85.1 cm)
Overall Depth:	72.2" (183.4 cm)
Overall Height (down):	43.7" (111 cm)
Overall Height (raised):	57.7" (146.6 cm)
Max. Vertical Lift:	n/a
Box 1	80cm x 86cm x 83cm
Box 2	n/a
Net weight	53 kg
Gross Weight	59 kg