

The Hoyer® Presence® lift has been engineered and designed to handle almost any resident-handling task. Along with the Stature lift, the Presence can lift residents from the floor as well accommodate transfers to higher surfaces. The 500 lbs safe working load and greater spatial area make this lift a perfect fit for larger residents. The Presence comes standard with a powered base and is available with a scale option.





## Smart Monitor

The unique Smart Monitor control box is programmable and will accurately record and display vital lift utilization and service data. An easy to read LCD display provides the valuable information at the touch of a button.



## Interchangeable Cradle Bars

Part of the restriction of using certain lifts are the slings that you can use on the product. The Presence however, gives you a choice of sling systems.

| SPECIFICATION                           | IMPERIAL          | METRIC  |
|---|-------------------|---------|
| Safe Working Load                       | 500 lbs           | 227 kgs |
| Maximum Overall Length                  | 57.5"             | 1460 mm |
| Minimum Overall Length                  | 55.1"             | 1400 mm |
| Maximum Overall Height                  | 81.3"             | 2065 mm |
| Minimum Overall Height                  | 57.5"             | 1460 mm |
| Spreader Bar Maximum Height             | 75.2"             | 1910 mm |
| Spreader Bar Minimum Height             | 16.9"             | 430 mm  |
| Height at Maximum Reach                 | 51.5"             | 1310 mm |
| Reach at Maximum Height                 | 25.6"             | 650 mm  |
| Reach at Minimum Height                 | 15.7"             | 400 mm  |
| Maximum Reach                           | 36.6"             | 930 mm  |
| Turning Radius                          | 64.2"             | 1630 mm |
| Legs Open - Internal Width              | 40.1"             | 1020 mm |
| Legs Closed - External Width            | 26.4"             | 670 mm  |
| Front Twin Casters                      | 4.0"              | 100 mm  |
| Rear Braked Casters                     | 4.0"              | 100 mm  |
| Overall Height of Legs                  | 4.7"              | 120 mm  |
| With Low Leg Option                     | 2.4"              | 60 mm   |
| Ground Clearance                        | 1.4"              | 35 mm   |
| * Deach contar of an used or har to the | front of the most |         |

\* Reach = center of spreader bar to the front of the mast

#### 4 POINT POSITION CRADLE (measurement to top of location pins)

| Cradle Max. Height          | 61.0" | 1550 mm |
|-----------------------------|-------|---------|
| Cradle Min. Height (usable) | 15.4" | 390 mm  |

| WEIGHT | IMPERIAL | METRIC |
|--------|----------|--------|
| Total  | 88.0 lbs | 40 kgs |

## ELECTRIC SHOCK PROTECTION

Charger – Class II Lift - Internal power source

#### DEGREE OF SHOCK PROTECTION

Charger – Type B Lift – Type B Intended operating environment: >+5°C <+40°C Outside this environment functionality and safety may be compromised

## **ELECTRICAL SPECIFICATIONS**

Battery – 24vDC rechargeable sealed lead acid type 3.2A hours Charger rated input - 100-240V AC 24 VDC 50/60Hz Charger rated output – 29.5 VDC, Max. 19 W



Hoyer Presence lift is designed to fit under beds, stretchers, and other furnishings with low



# Low Leg Option The Low Leg Option Kit for the

clearance. Total leg height is 2.4". Part #: 39001310

# Foot Pad

Ergonomic foot pad to assist with initiating movement.