





The Essential Collection from Pride®
Power Lift Recliners provides exceptional
value and outstanding comfort.
Several fabric options will match any
décor while delivering rest and worryfree operation.













	Medium	Medium	Medium
SPECIFICATIONS	LC-250	L-158	LC-105
	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'2"- 5'10"	3-Position Chaise Lounger 5'2" - 5'10"
Weight Capacity	375 lbs.	325 lbs.	325 lbs.
Seat-to-Floor	20"	18"	18"
Seat Depth	20"	18"	19"
Seat Width	22"	21"	20"
Top of Back to Seat	26"	25"	28"
Back Style	Blown Fiber Split Back	Blown Fiber Button Back	Sewn Pillow Back
Fully Padded Chaise	Yes	Yes	Yes
Heat & Massage	No	No	No
Seating Type	Total Comfort	Sinuous Spring	Sinuous Spring
Distance From Wall	20"	18"	17"
Overall Chair Width	34"	33"	33"
Chair Weight	118 lbs.	102 lbs.	96 lbs.
Warranty	2-Year Limited	1-Year Limited	1-Year Limited

LEXIS STA-KLEEN® BENEFITS

- Easy-to-use hand control adjusts your chair to find your perfect position
- Large dual pockets offer ample storage space that is easily accessible
- Smooth and durable lift and scissor mechanism for long-lasting and quiet performance
- Removable back lets you adjust chair padding as needed
- Integrated battery backup ensures chair is operational when the power goes out1
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Lifetime warranty on steel lift frame
- Fire retardant (California specifications; bulletin 117-2013)

OPTIONS*

- · Residential cleaning kit
- Protective pad
- · Non-skid leg levelers
- *Additional charges apply for all options.
- 1. Batteries not included



L-250 FABRIC SELECTION















OPTIONAL

LC-158 FABRIC SELECTION

EL CAMINO











LC-105 FABRIC SELECTION

MICRO-SUEDE











OPTIONAL

LEXIS URETHANE®



Look for the Quick Ship Indicator for the shortest lead times available



3-POSITION FULL RECLINE The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.