

<b>Model</b>	Go-Go Elite Traveller®
<b>Part Number</b>	3-wheel: SC40E 4-wheel: SC44E
<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b> <sup>1</sup>	Up to 4 mph
<b>Ground Clearance</b> <sup>3</sup>	1.38" at motor / 0.75" at anti-tips
<b>Turning Radius</b> <sup>3</sup>	3-wheel: 33" 4-wheel: 44"
<b>Overall Length</b> <sup>3</sup>	3-wheel: 37" 4-wheel: 39.5"
<b>Overall Width</b> <sup>3</sup>	19.5"
<b>Front Tires</b>	3-wheel: 2" x 8" solid 4-wheel: 2" x 7" solid
<b>Rear Tires</b>	2.5" x 8" solid
<b>Suspension</b>	Limited
<b>Range Per Charge</b> <sup>1,2</sup>	3W: Up to 8.5 miles (12 AH) 4W: Up to 8.2 miles (12 AH) 3W: Up to 12.4 miles (18 AH) 4W: Up to 12 miles (18 AH)
<b>Total Weight Without Batteries</b>	3-wheel: 76 lbs. 4-wheel: 80 lbs.
<b>Heaviest Piece When Disassembled</b> <sup>3</sup>	3-wheel: 29 lbs. (front section) 4-wheel: 34 lbs. (front section)
<b>Standard Seating</b>	Type: Foldable, compact, vinyl-covered Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17" Depth: 17"
<b>Motor Output</b>	24-volt DC motor, sealed mini-transaxle
<b>Battery Requirements</b> <sup>4</sup>	Size: (2) 12 AH Weight: 9 lbs. each, 20 lbs. per battery pack or Size: (2) 18 AH Weight: 14 lbs. each (optional), 29.5 lbs. per battery pack
<b>Battery Charger</b>	Off-board, 2 amp (12 AH), Off-board, 3.5 amp (18 AH)

### HCPCS Code: K0800 (Group 1, Standard)

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> Battery weight may vary based on battery manufacturer.

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GOGOET\_061819

**GOGO**  
ELITE TRAVELLER<sup>®</sup>

**THE  
QUICKEST  
WAY TO  
TRAVEL**



**LIVE YOUR BEST<sup>®</sup>**

pridemobility.com |    

# THE QUICKEST WAY TO TRAVEL



Two sets of easily changeable, red and blue colored panels.



Feather-touch disassembly for quick transport or storage.



Charger port located conveniently on the tiller.

The compact design of the Go-Go Elite Traveller® allows you to easily maneuver in tight spaces while providing stable outdoor performance. With a 300 lb. weight capacity, a charger port in the tiller, one-hand feather-touch disassembly and maximum speeds up to 4 mph, it's easy to go where you want to go!

## Additional Features

- Delta tiller with charger port to conveniently charge your scooter
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces
- Front frame-mounted seat post offers maximum stability
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires
- Removable footboard basket for additional storage (4-wheel version only)

## Colors (both color panel sets included)



Red



Blue

## Accessories\*

- Cane clip attachment
- Cane/crutch holder (single, double)
- Cup holder
- Forearm crutch holder
- 50" Lap belt
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- extra 12AH battery box
- extra 18AH battery box

## Military Patches



\* Only one rear accessory can be used at a time.