

# Jazzy<sup>®</sup>

## ELITE 14

The Jazzy<sup>®</sup> Elite 14 is designed for excellent indoor maneuverability and outstanding outdoor capability. The Jazzy Elite 14 features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns around corners. Large front wheels offer extra absorption over obstacles for superb climbing capabilities along with front anti-tip wheels for transitioning over grade changes and handling various terrains with ease.



### LARGE 14" DRIVE WHEELS



### PERFORMANCE-MATCHED ELECTRONICS



### COLOR-THROUGH SHROUD MEANS SCRATCHES ARE LESS VISIBLE

## AMERICA'S FAVORITE POWER CHAIR

Learn more online.



[gojazzy.com](http://gojazzy.com)

# Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	3"
<b>Turning Radius</b>	24"
<b>Overall Length</b>	35.5"
<b>Overall Width</b>	24.25"
<b>Drive Wheels</b>	14" solid
<b>Front Wheels</b>	3" solid anti-tips
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, front-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	50A, Dynamic Shark controller
<b>Battery Charger</b>	Off-board, 5 amp
<b>Per-Charge Range</b>	Up to 9.8 miles
<b>Battery Requirements</b>	(2) 12 volt, NF-22
<b>Battery Weight</b>	37 lbs. each
<b>Base Weight</b>	99 lbs.
<b>Standard Seat Weight</b>	52 lbs. (high-back)
<b>Standard Seat</b>	18" x 18-20", 20" x 18-20" and 22" x 20-22"
<b>Max. Seat Dimensions</b>	22" x 22"
<b>Warranty</b>	5-year limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

## Colors



Jazzy Red



Jazzy Blue

## Models

<b>Jazzy Elite 14-2S-SS</b>	Code: K0822
<b>Jazzy Elite 14-2S-C</b>	Code: K0823

\*Additional charges may apply.

# Features and Benefits

- Color-through shroud
- Depth adjustable, limited recline comfort high back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## Seating Options\*

- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16"x16" or 20"x20"
- Synergy® Seat

## Options

- Bell shaped foot platform
- Elevating leg rests
- Heel loops
- Height adjustable foot platform bracket
- Residual limb support
- Swing- away legrests
- Swing-away joystick

## Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

## Recommended Lifts

- Backpacker® Series
- Outlander Full Platform
- Commander 450



Find a compatible lift for your vehicle online:  
[PRIDEMOBILITY.COM/LIFT-FINDER](http://PRIDEMOBILITY.COM/LIFT-FINDER)



(US) 800-800-8586 • (Canada) 888-570-1113 • [pridemobility.com](http://pridemobility.com) •



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.