

SATURN LIFT CHAIR



Saturn CH4005 (Blue) CH4006 (Brown)

The Saturn lift chair is an infinite position lift chair available with brown or blue upholstery. While it comes fitted with rubber pads to prevent floor damage and provide stable positioning, the 4 supplied locking casters can be installed to ease repositioning the chair. Additionally, the Saturn comes equipped with an adjustable tray and a cup holder for user convenience.

With a Wide/Average size rating, this chair will accommodate people of wider widths and average height.

- Infinite position recline
- Four locking swivel casters included
- Adjustable tray and cup holder



The spring assisted footrest on all models will flip up if it comes in contact with an obstruction when lowering to minimize injury or damage.



Supplied locking swivel casters can be installed on all models, except Earth and Mars which have 3" casters. Ideal for users who may want to move their chair occasionally or rental units.



Our Saturn and Mercury models are upgraded to include a cup holder and adjustable tray for the convenience of the user. Mercury also includes leather upholstery and heat & massage.

FEATURES:

- Lounger style that fully supports the legs (no gap between footrest & seat)
- Removable back for easy handling
- Easy to use hand control
- UL / CUL Approved
- Maximum seat elevation with superb stability in full lift position
- Wide construction provides greater side to side stability
- Top steel base for durability & cross supports
- Engineered bushings & washers at pivot points provide quiet smooth operation & durability
- Direct drive activation synchronizes right & left sides for easy opening
- Low voltage motor, 24 volt
- Battery Back Up



MODEL NAMES	SATURN
Model numbers:	CH4005 (Blue) CH4006 (Brown)
Weight Capacity:	352 lbs (160 kg)
Seat Width:	21.1" (53.6 cm)
Seat Depth:	19.5" (49.5 cm)
Back Height:	28.5" (72.4 cm)
Seat Height:	19.5" (49.5 cm)
Overall Width:	31.2" (79.2 cm)
Overall Depth:	73.3" (186.2 cm)
Overall Height (down):	43.3" (110 cm)
Overall Height (raised):	58.5" (148.6 cm)
Max. Vertical Lift:	n/a
Box 1	81cm x 83cm x 79cm
Box 2	n/a
Net weight	54 kg
Gross Weight	60 kg